

8000GC/MSD Gas chromatography Mass spectrometry



NWspec MS Factory

Zhanye Road 18th, SIP, Suzhou, Jiangsu, China, 215000

www.nwspec.com






CORPORATION

PROFILE

NWSpec is a leading high-tech enterprise in China engaged in the research and development, production, sales, and service of chemical analysis and medical detection mass spectrometer. We focus on providing users with fast, easy-to-use, and highly robust laboratory mass spectrometers and portable mass spectrometers for various application scenarios, and have accumulated advanced technologies at an international level and multiple core patents.

Existing products: LC-TQMS, GC-TQMS, GC-MS, and Portable GC-MS.



HARDWARE

features

NWSpec 8000 series products is a high-end GC-MS independently developed based on our own patented technology, with a detection limit of better than 10^{-14} g, which is at the top level of similar products in the world. It can be widely used in high demand fields such as scientific research, pesticide residue detection, environmental monitoring, and metabolomics research.

All metal conjugate quadrupole

The cold rod design does not require heating and lifelong maintenance. It adopts the most advanced segmented design, with detachable pre-quadrupoles that can automatically adjust DC voltage, eliminate edge field effects, prevent ion leakage, and thus achieve higher ion transmission efficiency.



High temperature inert ceramic ion source

Efficient ionization, reducing pollution, equipped with two long lifespan special material filaments, providing double usage time, and high-precision temperature control for all lenses. Easy to disassemble and clean.



Independent pre-quadrupole

Applying RF voltage, combined with pre-quadrupole technology, effectively eliminates neutral fragment interference and reduces noise; Simultaneously optimize the ion source and quadrupole transition electric field.




Discontinuous discrete electron multiplier

Longer service life, stronger resistance to pollution and water.



Equipped with a dual-cavity ultra-high vacuum system as standard



To achieve a better vacuum environment, a dual-cavity vacuum chamber has been developed to enable individual exhaust for ion source and quadrupole part. A ultra-high vacuum environment has the following advantages:



Providing the whole MS system with a high vacuum level can not only reduce the contamination of the ion source, but also maintenance operation.




Fast turn-on enabled. Users can start their experiment within 30min after startup.



Maximum column flow supports up to 5mL/min (10mL/min available on flagship model). More advanced applications can be implemented by using larger-diameter GC columns.



Electron multiplier in high vacuum status benefits ultra-high sensitivity test. GC direct injection, OFN method IDL is better than 10^{-14} g.



定性表 - 方法: -								
	物质名称	CAS	保留时间	定量离子	参考离子	实际响应比%	参考响应比%	响应比偏差
1	特丁硫醇		16.464	231.0→128.9	231.0→174.9, 231.0→202.9	100.0, 71.6, 56.7	100.0, 71.6, 56.7	30, 30, 30
2	地虫硫醇		16.542	137.1→109.1	246.0→137.1, 246.0→109.1	100.0, 65.3, 91.1	100.0, 65.3, 91.1	30, 30, 30

Through the qualitative table function, the rapid qualitative of multi-substance and multi-data is realized

数据结果统计											
物质名称	序号	文件名	Rt	Rt 平均值	Rt KSD [%]	响应值	响应平均值	响应值 KSD [%]	含量	含量平均值	含量 KSD [%]
六六六	1	痕量性有机氯100ppb_6(1).nxd	8.855	8.848	0.038	1063.28	1030.32	3.02	102.67	99.72	2.79
	2	痕量性有机氯100ppb_6(2).nxd	8.848			1067.78			103.07		
	3	痕量性有机氯100ppb_6(3).nxd	8.847			997.81			96.81		
	4	痕量性有机氯100ppb_6(4).nxd	8.847			1040			100.58		
	5	痕量性有机氯100ppb_6(5).nxd	8.848			1001.2			97.11		
	6	痕量性有机氯100ppb_6(6).nxd	8.844			1011.85			98.06		
六氯苯	1	痕量性有机氯100ppb_6(1).nxd	9.025	9.018	0.042	749.48	726.05	2.64	113.17	109.58	2.68
	2	痕量性有机氯100ppb_6(2).nxd	9.019			712.16			107.46		
	3	痕量性有机氯100ppb_6(3).nxd	9.018			748.09			112.96		
	4	痕量性有机氯100ppb_6(4).nxd	9.015			707.89			106.8		
	5	痕量性有机氯100ppb_6(5).nxd	9.018			729.21			110.07		
	6	痕量性有机氯100ppb_6(6).nxd	9.015			709.46			107.04		
六六六	1	痕量性有机氯100ppb_6(1).nxd	9.307	9.299	0.046	725.06	684.81	3.88	104.67	99.71	3.28
	2	痕量性有机氯100ppb_6(2).nxd	9.297			703.5			102.01		
	3	痕量性有机氯100ppb_6(3).nxd	9.297			659.7			96.62		
	4	痕量性有机氯100ppb_6(4).nxd	9.297			662.17			96.83		
	5	痕量性有机氯100ppb_6(5).nxd	9.297			692.54			100.67		
	6	痕量性有机氯100ppb_6(6).nxd	9.297			665.91			97.39		
六六六	1	痕量性有机氯100ppb_6(1).nxd	9.457	9.45	0.04	749	706.95	4.31	104.41	99.45	3.61
	2	痕量性有机氯100ppb_6(2).nxd	9.447			684.49			97.99		
	3	痕量性有机氯100ppb_6(3).nxd	9.45			683.23			97.64		
	4	痕量性有机氯100ppb_6(4).nxd	9.447			737.92			103.11		
	5	痕量性有机氯100ppb_6(5).nxd	9.451			698.74			98.49		
	6	痕量性有机氯100ppb_6(6).nxd	9.451			668.3			94.9		
六六六	1	痕量性有机氯100ppb_6(1).nxd	9.89	9.88	0.048	686.31	635.3	6.01	103.41	96.9	5.02
	2	痕量性有机氯100ppb_6(2).nxd	9.876			663.54			100.5		
	3	痕量性有机氯100ppb_6(3).nxd	9.879			615.93			94.43		
	4	痕量性有机氯100ppb_6(4).nxd	9.879			654.53			99.35		
	5	痕量性有机氯100ppb_6(5).nxd	9.879			597.17			92.04		
	6	痕量性有机氯100ppb_6(6).nxd	9.879			584.31			91.67		
七氯	1	痕量性有机氯100ppb_6(1).nxd	10.993	10.987	0.031	151.06	129.74	9.41	120.45	109.46	5.74
	2	痕量性有机氯100ppb_6(2).nxd	10.987			132.39			110.83		
	3	痕量性有机氯100ppb_6(3).nxd	10.986			121.11			106.02		
	4	痕量性有机氯100ppb_6(4).nxd	10.987			132.23			110.75		
	5	痕量性有机氯100ppb_6(5).nxd	10.987			125.48			107.27		
	6	痕量性有机氯100ppb_6(6).nxd	10.983			116.16			102.47		

It has the function of data aggregation and statistics




Soil and sediment Determination of semi-volatile organic compounds

Summary of methodology

Refer to HJ 834-2017 Determination of semi-volatile organic compounds in soil and sediment - Gas chromatography-mass spectrometry

Test Method

Injection Volume:	No shunts
Inlet:	1.0μL
Column(s):	280°C
Column Flow:	DB-5MS 30m×0.25mm×0.25μm constant current mode
Column box heating procedure:	1ml/min
Ion source:	35°C (2min) 15°C/min 150°C (5min) 3°C/min 290°C (2min)
Transmission Lines:	280°C
Solvent Delay:	280°C
Scan Mode:	3min
SCAN Mode Scan Range:	Full Scan (SCAN) and Selective Ion (SIM) 35~450amu



Total ion chromatogram

Ambient air Determination of volatile organic compounds

Summary of methodology

Refer to HJ 644-2013 Ambient air - Determination of volatile organic compounds - Adsorption tube sampling-thermal desorption/gas chromatography-mass spectrometry

Thermal desorption conditions


Temperature setting: adsorption tube 300 °C; Transmission line 200°C; Injection valve 160°C; Cold trap -20 °C
Focus tube: 300°C; The aging temperature is 310°C
Time Setting: Purging time 251s
Injection time: 200s
Aging time: 800s

Gas phase conditions

Injection Mode: shunt
Shunt ratio: 10:1
Inlet: 220°C
Carrier gas: Helium
Column(s): DB-624 60 m×0.25 mm×1.4 μm constant current mode
Column Flow: 1ml/min
Column box heating procedure: 40°C (3 min) 10°C/min 200°C (6 min)

Mass spectrometry conditions

Ion source: 240°C
Transmission Lines: 250°C
Solvent Delay: 7.5 min
Scan Mode: Full Scan (SCAN) and Selective Ion (SIM)
SCAN Mode Scan Range: 35~300amu



Total ion current plot of volatile organic compounds


Determination of volatile organic compounds in water quality

Summary of methodology

Refer to HJ 639-2012 Determination of Volatile Organic Compounds in Water Quality Purging and Trapping/Gas Chromatography-Mass Spectrometry

Test Method

Purging time: 11 min
Desorption temperature: 250°C
Desorption time: 1 min
Injection Mode: Split injection
Shunt ratio: 20:1
Inlet Temperature: 200°C
Column(s): DB-624 30m×250μm×1.4μm constant current mode
Column Flow: 1ml/min
Column box heating procedure: 38°C (1.8min) 10°C/min 120°C 15°C/min 240°C (5 min)
Ion source: 300°C
Transmission Lines: 250°C
Solvent Delay: 1.5 min
Scan Mode: Full Scan (SCAN) or Selective Ion (SIM)
Scan Range: 35~270amu



20 µg/L standard solution selective ion chromatogram

Drinking water testing


Determination of geosmin and 2-methylisoquinol in drinking water

Summary of methodology

The new sanitary standard for drinking water (GB5749-2022) has added two detection indicators of geosmin (GSM) and 2-methylisocamphenol (2-MIB). GB/T5750.8-2023 specifies the use of headspace solid-phase microextraction gas chromatography mass spectrometry for detection. In this report, the Anyipu 7700 GC/MS instrument is used. The method development of geosmin and 2-methylisoquinol in drinking water was carried out by headspace solid-phase microextraction gas chromatography-mass spectrometry, and the spirit of the instrument was evaluated. Sensitivity, stability, accuracy.

Mass spectrometry instrument conditions

Inlet Temperature:	250°C
Injection Method:	No shunts
Carrier gas:	helium
Flow rate:	1.00ml/min
Program Heating:	Keep at 60°C for 2.5min, raise the temperature at 8°C/min to 250°C for 5min
Ion source temperature:	250°C
Transmission Line Temperature:	280°C
Scan Mode:	Choosing an Ion Scan (SIM)
2-Methylisocamphenol characteristic ions:	95 (quantitative), 107 (qualitative 1), 135 (qualitative 2)
Geosmin characteristic ions	112 (quantitative), 125 (qualitative)




Analysis of orthosilicate components

Methyl orthosilicate sample method

Columns	DB-5MS 30m×0.25mm×0.25μm
Heating procedure:	45°C (1min) 20°C/min 250°C (5min)
Inlet Mode:	shunt
Shunt ratio:	10: 1
Inlet Temperature:	250°C
Column Flow Rate:	1.0mL/min (constant flow)
Injection Volume:	1μL
Ion source temperature:	240°C
Transmission Lines:	250°C
Solvent Delay:	0.1min
Scan Mode:	Scan
Time Segmentation:	0.1~2.5min 15~400amu, 2.5~3.8min off multiplier, Filament, 3.8~16.25min, 35~400amu

Ethyl orthosilicate sample method 0.1~4.3min 15~400amu, 4.3~5.9min off multiplier,
Time Segmentation: Filament, 5.9~16.25min, 35~400amu,
Other gas chromatography mass spectrometry conditions are the same as
methyl orthosilicate methods



Total ion current plot of the sample impurities

Food testing


Liquor composition analysis

Gas phase conditions

Injection Method:	Direct injection
Injection Volume:	1μL
Inlet Mode:	shunt
Shunt ratio:	10: 1
Inlet Temperature:	250°C
Column Flow Rate:	1.0mL/min (constant flow)
Column(s):	DM-WAX 60m×0.25mm×0.25μm
Column box heating procedure:	37°C (2min) 3°C/min 70°C (1min) 6°C/min 130°C (0min) 10°C/min 220°C (10min)

Mass spectrometry conditions

Ion source temperature:	240°C
Transmission Lines:	240°C
Solvent Delay:	0.1min
Filament off time:	8.3min~9.5min
Scan Mode:	Scan
Scan Range:	25~200amu
Ionization Energy:	70ev
Filament Current:	100μA
EM Voltage:	900V



Liquor total ion current (TIC) diagram

Chemical product testing


Determination of qualitative amines

Summary of methodology

The GCMS8000 gas chromatography-mass spectrometer was used to characterize the given samples. The qualitative method of the sample was finalized. In the Scan scanning method, the mass-spectra obtained by The components are: methanol (CAS number 67-56-1), ethylenediamine (CAS number 107-15-3), ethanolamine 10-85-0), diethylenetriamine (CAS 111-40-0), hydroxyethylethylenediamine (CAS 111-41-1), N-aminoethylpiperazine -31-8) and hydroxyethylpiperazine (CAS number: 103-76-4)

Test Method

Injection Mode:	shunt
Shunt ratio:	200: 1
Inlet Temperature:	280°C
Injection volume	0.2μL
Column(s):	DB-17MS 30m×0.25mm×0.25μm constant current mode
Column Flow:	1ml/min
Column box heating procedure:	60°C (1min) 30°C/min 250°C (2min)
Ion source:	240°C
Transmission Lines:	280°C
Solvent Delay:	0.1min
Scan Mode:	Full Scan (SCAN)
Scan Range:	30~300amu



Total ion current plot of amine sample full sweep (Scan).

Electrical and electronic product testing

Determination of 4 phthalates

Summary of methodology


Refer to GB/T 29786-2013 Determination of dimethyl phthalate in electrical and electronic products - Gas chromatography mass spectrometry

Chromatographic parameters

Injection Mode:	No shunts
Purging time:	1.0min
Purge Flow:	50ml/min
Injection Volume:	1.0μL
Inlet Temperature:	280°C
Column(s):	HP-5MS 30m×0.25mm×0.25μm constant current mode
Column Flow:	1ml/min
Column box heating procedure:	60°C (1min) 30°C/min 280°C (6min)

Mass spectrometry parameters

Ion source:	300°C
Transmission Lines:	290°C
Ionization Energy:	70eV
Filament Current:	200μA
EM voltage:	990V
Solvent delay	6min
Scan Mode:	Full Scan (SCAN) and Selective Ion (SIM)
SCAN Mode Scan Range:	45~300amu



Concentration 5.0 μg/ml 4 phthalates full sweep total ion current plot


Polybrominated biphenyls and polybrominated diphenyl ethers GCMS schemes

Summary of methodology

Refer to GB/T 26125-2011 Determination of six restricted substances (lead, mercury, cadmium, hexavalent chromium, polybrominated biphenyls and polybrominated diphenyl ethers) in electrical and electronic products

Test Method

Injection Mode:	Pulsed splitless injection
Pulse Pressure:	20psi
Pulse time:	0.75min
Purge flow	60ml/min
Purging time:	1min
Injection Volume:	1μL
Column(s):	DB-5HT, 15m×0.25mm×0.1μm
Column Flow:	3ml/min
Inlet:	280°C
Ion source:	350°C
Transmission Lines:	340°C
Solvent Delay:	2.5min
Multiplier voltage:	Tuning voltage +0.2KV
Scan Mode:	Full Scan (SCAN) or Selective Ion (SIM)
Scan Range:	100~1000amu



Chromatograms of 17 PBBs and PBDEs mixed with standard 5 mg/L



Peking university




A Center for Disease Control in Guangdong



A large third party in Guangxi



Nanjing Agricultural University



A environmental monitoring station
in Inner Mongolia



An agricultural product inspection center
in Qingdao



Zhejiang Electronic Technology Co., Ltd



Zhejiang Research Institute of Chemical Industry