



NWSPEC (Xiongan) Co., Ltd.  
Add. 15-051# Xiongan Pilot Free  
Trade Zone Area, Hebei, China

## ICP Technical Specifications

## Table of Contents

1.	Overview .....	3
1.1.	About Us .....	3
1.2.	Main structure of ICP spectrometer .....	3
1.3.	Brief description of the working principle of ICP spectrometer .....	4
1.4.	Range of applications .....	4
1.5.	Main parameters .....	5
1.6.	Main technical indicators .....	6
1.7.	Product standards .....	6
2.	Working Environment & Facility Requirements .....	7
2.1.	Room Size: The effective usable area of the room should be at least 15 square .....	7
2.2.	Placement: The rear side of ICP10 spectrometer should be at least 1000 .....	7
2.3.	Supporting Facilities: Reagent cabinet or reagent rack, gas cylinder cabinet, desks .....	7
2.4.	Operating Environment .....	7
2.5.	Dust Control Measures .....	7
2.6.	Stable power supply system .....	7
2.7.	Power Requirements .....	8
2.8.	Grounding Requirement .....	9
2.9.	Gas Requirements .....	9
2.10.	Cooling Water Requirements .....	9
2.11.	Exhaust Requirements .....	9
3.	Installation .....	10
3.1.	Unpacking and Inspection .....	10
3.2.	Installation .....	10
3.3.	Acceptance .....	11
4.	Hardware Operation .....	12
4.1.	Startup Procedure .....	12
4.2.	Ignition Procedure .....	12
4.3.	Flame shut-off Procedure .....	12
4.4.	Shutdown Procedure .....	12
5.	Software operation Procedure .....	13
5.1.	Zero Scan .....	13
5.2.	New Method .....	13
5.3.	Select Method .....	13
5.4.	Auto Search .....	13
5.5.	Auto Attenuate .....	13
5.6.	Standard .....	14
5.7.	Sample .....	14
5.8.	View Test Results .....	14
	15	
5.6.	Standard .....	15
5.7.	Sample .....	15
5.8.	View Test Results .....	15
6.	Common Fault Resolution .....	15-18

## 1. Overview

### 1.1. About Us

NWspec self-developed ICP spectrometer research and development, production, sales, and technical services, Inductively Coupled Plasma Optical Emission Spectrometer (ICP-OES). This instrument builds upon traditional ICP-OES technology while incorporating the latest international advancements, achieving full automation, high precision, high sensitivity, and exceptional stability.

With a dedicated R&D team committed to innovation and collaboration, the company continuously pushes technological boundaries. Backed by extensive industry experience and solid expertise, it ensures high-quality and professional service.

ICP10 & ICPO3 is an advanced analytical instrument that integrates optical, mechanical, electrical, computing, and analytical technologies. It features high-speed testing, a wide measurement range, and highly accurate and reliable analytical results. Leveraging computer technology, the instrument offers intelligent operation, graphical and textual data display, and efficient data acquisition and processing. These capabilities place it at the forefront of domestic technology, making it an ideal analytical tool for various industries.

### 1.2. Main structure of ICP spectrometer

#### 1.2.1. Sample Introduction System

The sample introduction system is responsible for delivering the sample into the plasma for analysis. It consists of a nebulizer, a spray chamber, and a carrier gas system. The liquid sample is nebulized into fine droplets and transported to the plasma torch by the carrier gas (argon) for excitation.

#### 1.2.2. Excitation Source

The sample introduction system is responsible for delivering the sample into the plasma for analysis. It consists of a nebulizer, a spray chamber, and a carrier gas system. The liquid sample is nebulized into fine droplets and transported to the plasma torch by the carrier gas (argon) for excitation.

#### 1.2.3. Optical Dispersion System

The main function of this system is to disperse and separate the multi-element spectral emissions from the plasma. A grating spectrometer is used to separate different wavelengths of light and focus them onto the detector. The optical system includes collimating mirrors, diffraction gratings, and dispersion prisms, ensuring high resolution and low stray light.

1.2.4. Detection and Control System

The main function of this system is to disperse and separate the multi-element spectral emissions from the plasma. A grating spectrometer is used to separate different wavelengths of light and focus them onto the detector. The optical system includes collimating mirrors, diffraction gratings, and dispersion prisms, ensuring high resolution and low stray light.

1.2.5. Output System

The output system includes computer software and data interfaces for displaying results, storing data, and generating reports. Users can configure operating parameters, perform spectral analysis, conduct quantitative calculations, and export experimental data through the software interface to meet various analytical requirements.

1.3. Brief description of the working principle of ICP spectrometer

The ICP spectrometer excites sample elements using high-temperature plasma, causing them to emit characteristic spectra for qualitative and quantitative analysis. The sample is atomized by the sample introduction system, mixed with carrier gas, and introduced into the plasma torch. A high-frequency electromagnetic field generates a high-temperature plasma, which evaporates, dissociates, and excites the sample to a high-energy state.

As the excited atoms and ions return to their ground state, they release light at characteristic wavelengths. This light is dispersed by the optical system, detected, and converted into electrical signals, which are then processed and analyzed. Finally, the analysis results are displayed, stored, and reported using computer software.

1.4. Range of applications

- 1.4.1. **Metallurgical Industry:** Steel and its alloys, including carbon steel, cast iron, alloy steel, high-purity steel, silicon and iron alloys, ferroalloys, and rare earth ferroalloys.

- 1.4.2. **Non-Ferrous Metals and Their Alloys:** Including non-ferrous metals (high-purity metals) and their alloys, rare metals and their alloys, precious metals, rare earth elements and their alloys, as well as their compounds.
- 1.4.3. **Water Quality Samples:** Drinking water, surface water, wastewater, mineralized water, and electroplating solutions and wastewater.
- 1.4.4. **Environmental Samples:** Soil, atmospheric dust, fly ash, and solid waste.
- 1.4.5. **Geological and Mining Samples:** Geological samples, ores, and minerals.
- 1.4.6. **Inorganic Non-Metallic Material Analysis.**
- 1.4.7. **Medical and Health:** Detection of Ca, Mg, Cu, Pb, Cd, Fe, Mn, Se, and Zn requires only 20  $\mu\text{L}$  of blood or 0.1000 g of hair.
- 1.4.8. **Chemical and Industrial Products:** Chemical reagents (hydrochloric acid, nitric acid, sulfuric acid, phosphoric acid, hydrofluoric acid, and other synthetic reagents), industrial chemicals (e.g., saturated sodium chloride solution), inorganic materials, cosmetics, oils (gasoline, diesel, mineral oils, marine oils, petroleum coke, and raw materials), petroleum catalysts, urea solutions, and coatings.
- 1.4.9. **Agricultural and Food Products:** Analysis of metal elements in grains, oils, seafood, food, and beverages.
- 1.4.10. **Animal, Plant, and Biochemical Samples:** Plants, traditional Chinese medicine (TCM), animal tissues, and biochemical samples.
- 1.4.11. **Nuclear Industry Products:** Nuclear fuel analysis and nuclear materials.
- 1.4.12. **Battery Materials:** Lithium carbonate, lithium cobalt oxide (LCO), lithium manganese oxide (LMO), nickel-cobalt-manganese (NMC) ternary materials, and lithium iron phosphate (LFP).
- 1.5. Main parameters
  - 1.5.1. **High-quality optical components with an advanced control system:**  
Ensures precise scanning and positioning to achieve a high signal-to-noise ratio in test results.
  - 1.5.2. **Key components sourced from France, Germany, and Japan:** Guarantees high accuracy and sensitivity of the instrument.
  - 1.5.3. **Minimal matrix effect:** 99% of samples do not require matrix separation.
  - 1.5.4. **Wide measurement range:** Capable of analyzing from ultra-trace to major concentrations.
  - 1.5.5. **High measurement accuracy:** Relative standard deviation (RSD)  $\leq 1.5\%$ .
  - 1.5.6. **Fast processing speed:** Sequential scanning of more than 15 elements per minute.
  - 1.5.7. **Extensive element analysis:** Capable of detecting 70 metal elements and some non-metal elements.

- 1.5.8. **High stability:** Relative standard deviation (RSD)  $\leq 2.0\%$ .
- 1.5.9. **Low detection limits:** Meets national Class A standards.
- 1.5.10. **Supports both qualitative and quantitative analysis.**
- 1.5.11. **User-friendly analysis software:** Runs on **Windows 10 and Windows 11**, featuring a third-generation Chinese or English interface. The software is plug-and-play with no installation required, offering a powerful data processing system with multiple functions and flexible printing options.
- 1.6. Main technical indicators
  - 1.6.1. **Scanning range:**
    - **Grating 4320 lines/mm:** 180–442.5 nm
    - **Grating 3600 lines/mm:** 180–530.0 nm
    - **Grating 2400 lines/mm:** 180–800.0 nm
  - 1.6.2. **Wavelength indication error:**  $\leq \pm 0.03$  nm, **repeatability:**  $\leq 0.003$  nm
  - 1.6.3. **Minimum spectral bandwidth:** Mn 257.610 nm, better than national Class A standard,  $\leq 0.015$  nm
  - 1.6.4. **Detection limits (mg/L):**
    - **Zn 213.856 nm  $\leq 0.003$**
    - **Mn 257.610 nm  $\leq 0.002$**
    - **Ba 455.403 nm  $\leq 0.001$**
    - **Ni 231.604 nm  $\leq 0.01$**
    - **Cr 267.716 nm  $\leq 0.007$**
    - **Cu 324.754 nm  $\leq 0.007$**
  - 1.6.5. **Stability:** Relative standard deviation (RSD)  $\leq 2.0\%$
- 1.7. Product standards
 

ICP10 Inductively Coupled Plasma Emission Spectrometer, implementing the corporate standard Q/DXHKT0001-2017 of NWSPEC Co.,Ltd.

## 2. Working Environment & Facility Requirements

- 2.1. Room Size: The effective usable area of the room should be **at least** 15 square meters.
- 2.2. Placement: The rear side of ICP10 spectrometer should be **at least** 1000 mm away from the wall, and the right side should be **at least** 400 mm away from the wall.
- 2.3. Supporting Facilities: Reagent cabinet or reagent rack, gas cylinder cabinet, desks and chairs, **1–2 HP air conditioners**.
- 2.4. Operating Environment
  - **Ambient Temperature:** ( $25^{\circ}\text{C} \pm 3^{\circ}\text{C}$ ) with a temperature fluctuation rate within  $2^{\circ}\text{C}/\text{h}$ .
  - **Relative Humidity:**  $\leq 70\%$ , ideally **45-60%**. (Dehumidifiers should be used in high-humidity areas).
  - **Exhaust Requirements:** Doors and windows must be tightly sealed to prevent dust and acid mist intrusion.
  - **Other Requirements:** The room should be clean, dry, dust-free, free of corrosive gases, and free from vibrations.
  - **Grounding:** Proper grounding with a resistance of  $\leq 4\ \Omega$ .
- 2.5. Dust Control Measures

To ensure the long-term stable operation of the ICP spectrometer, the following dust prevention measures should be implemented to prevent dust ingress from affecting the performance of optical components, circuit systems, and other precision parts:

  - Wipe the instrument surface gently with a lint-free cloth moistened with a small amount of anhydrous ethanol or a neutral cleaner. Avoid using corrosive solvents.
  - Although the optical system is sealed, maintaining a clean laboratory environment is still essential.
  - The electronic components area should be protected from exposure to dusty environments. Periodic cleaning with low-pressure dry air is recommended to prevent dust accumulation.
  - Regularly inspect and clean the air vents to prevent dust blockages that could impact cooling performance.
- 2.6. Stable power supply system

To ensure the safe operation of the ICP instrument, the power supply line must have sufficient capacity. Otherwise, excessive voltage drop during

instrument operation may negatively impact its lifespan. As a precision measuring instrument, the ICP spectrometer requires a stable power supply, with voltage fluctuations generally not exceeding  $\pm 1\%$ . If the voltage variation surpasses this range, a precision AC voltage stabilizer with purification capabilities should be used. Standard voltage regulators often produce waveform clipping under high voltage conditions, leading to electrical impulses that can interfere with electronic computers, microprocessors, and phase-sensitive amplifiers, potentially causing malfunctions. A continuous sine wave power supply is essential to ensure the stable operation of these electronic circuits.

The power supply line for the instrument should ideally be connected directly to the distribution panel of the power transformer, separate from high-power equipment such as large motors, ventilators, air conditioners, and muffle furnaces. Sharing a power line with such equipment can lead to significant voltage fluctuations when these devices start, resulting in instrument instability. The grounding resistance for the spectrometer should be less than 4 ohms, while the computer grounding resistance should meet ASTM standards of less than 0.25 ohms to prevent electromagnetic interference.

During instrument operation, it is crucial to monitor power supply voltage fluctuations and avoid prolonged operation under overvoltage or undervoltage conditions. Overvoltage operation significantly reduces the lifespan of electron tubes to less than 20% of their normal lifespan. Conversely, undervoltage operation causes the filament temperature of electron tubes to drop excessively, leading to premature degradation of the electron-emitting material, further shortening the tube's lifespan. Additionally, large voltage fluctuations during instrument operation can cause instability in the high-frequency generator output power, significantly affecting measurement accuracy. Therefore, maintaining a high-quality power supply is essential for reliable instrument performance.

## 2.7. Power Requirements

ICP10 requires a single-phase AC 220V/50Hz power supply. If the voltage fluctuation exceeds  $\pm 1V$ , a 10KVA precision purified AC voltage stabilizer must be used (output accuracy  $\leq \pm 1\%$ , response time  $< 20-80ms$ ). The specific power requirements for each component are as follows:



- ICP10 Main Unit: Single-phase AC (220±1)V/50Hz, 2kVA (incoming wire cross-section ≥4mm<sup>2</sup>).
- Computer: Single-phase AC (220±1)V/50Hz, 500VA.
- Cooling Water Tank: Single-phase AC (220±1)V/50Hz, 1.5kVA.
- Exhaust System: Single-phase AC (220±1)V/50Hz, 20W (fan provided by the manufacturer).
- Circuit Breakers: Main switch: 60A , ICP: 32A , Cooling water tank: 20A , Axial fan: 10A

## 2.8. Grounding Requirement

The ICP main unit must be reliably grounded. The grounding substrate should meet the following minimum requirements:

- Copper Plate: At least 400mm (L) × 300mm (W) × 3mm (T)
- Copper Rod: At least  $\phi$  12mm × 1500mm
- Copper Tape: At least 20mm (W) × 0.5mm (T)
- The copper plate must be buried 1.5m underground.
- Pour 1kg of salt into the burial site to improve conductivity.
- Connect the copper tape from the buried copper plate to the instrument's indoor grounding point.
- The ground resistance must be less than 4Ω.
- When burying the copper plate, **remove large stones from the soil** to ensure proper contact with the earth.

## 2.9. Gas Requirements

Two bottles of argon gas with purity ≥99.99% (one bottle as a spare). The pressure reducing valve has an inlet pressure of 25 MPa and an outlet pressure adjustment range of 0–2.5 MPa. During use, set the outlet pressure to 0.2–0.25 MPa.

## 2.10. Cooling Water Requirements

When installing the smart cooling water tank, add 25 liters of distilled or purified water to the PVC tank of the cooling system.

## 2.11. Exhaust Requirements

The gas discharged from the plasma chamber contains argon, some metal vapors, and solvents. It is recommended to use  $\Phi$ 110 PVC pipe and  $\Phi$ 110 PVC adapters, along with an axial fan, to vent the exhaust outside. The outdoor section should be equipped with an elbow to prevent dust and rainwater from entering the duct.


### 3. Installation

#### 3.1. Unpacking and Inspection

After opening the wooden packaging box, remove the technical documents and check the instrument's components against the packing list to ensure all parts are present and undamaged. If any parts are missing or damaged, contact the manufacturer and the freight company immediately.

#### 3.2. Installation

- 3.2.1. Once all installation conditions are met, you can notify our technicians to come on-site for installation, commissioning, and training.
- 3.2.2. The ICP main unit is equipped with shock absorbers at the corners, which must be removed after unpacking. Once the instrument is positioned, engage the wheel brakes to stabilize it.
- 3.2.3. On the back of ICP10 spectrometer, from left to right, are the ICP main power cord, grounding terminal, argon gas inlet, and cooling water inlet and outlet.



- 3.2.4. ICP Mainframe Power Supply (① in the figure):  
Connect the ICP mainframe to a stabilized power supply, and then connect the stabilized power supply to a 32A circuit breaker. (If the voltage is stable, the stabilized power supply may not be necessary.) It is recommended to separate the phase wires supplying the ICP mainframe from those of the cooling circulation water tank and other ancillary equipment.
- 3.2.5. Ground Wire (② in the figure):

Use copper tape to connect the ICP mainframe to the outdoor grounding point.

3.2.6. Argon Gas Inlet (③ in the figure):

Use a  $\Phi 6$  PU pipe to connect the ICP mainframe to the argon gas cylinder.

3.2.7. Cooling Water Inlet and Outlet (④ in the figure):

Use a  $\Phi 10$  water pipe to connect the ICP mainframe to the water tank.

3.2.8. Auxiliary Power


Plug the auxiliary power supply into the plug board.

3.2.9. Communications Interface

Connect the ICP mainframe to the computer using a serial cable.

3.2.10. Installation of exhaust duct

Install the exhaust ducting according to the diagram below, connecting the axial fan as close to the outlet as possible.



3.3. Acceptance

After completing the installation of the instrument, the acceptance process can begin. Before proceeding, carefully read the instruction manual to understand the functions of all switches and buttons, and familiarize yourself with the instrument's operation methods and procedures before starting its operation.

## 4. Hardware Operation

### 4.1. Startup Procedure

- 4.1.1. Power on the cooling water tank
- 4.1.2. Turn on the axial fan
- 4.1.3. Turn on the argon
- 4.1.4. Turn on the ICP main power supply
- 4.1.5. Turn on the auxiliary control power cord
- 4.1.6. Turn on the main power on the right side of the instrument
- 4.1.7. Turn on the carrier gas switch, immerse the capillary tube in pure water or any solution, and check if the sample is fed smoothly
- 4.1.8. After completing the above steps, wait for the instrument to warm up for 3-5 minutes and then start the ignition

### 4.2. Ignition Procedure

- 4.2.1. Turn on the carrier gas for 10 seconds and then turn off the carrier gas
- 4.2.2. Turn on the plasma gas
- 4.2.3. Turn on the high voltage
- 4.2.4. Click the ignition button
- 4.2.5. Turn on the carrier gas

### 4.3. Flame shut-off Procedure

- 4.3.1. Shut down high pressure
- 4.3.2. Turn off the plasma gas
- 4.3.3. Shut off the carrier gas

### 4.4. Shutdown Procedure

- 4.4.1. Turn off the main power on the right side of the instrument
- 4.4.2. Turn off the auxiliary control power cord
- 4.4.3. Turn off the ICP main power supply
- 4.4.4. Turn off the argon
- 4.4.5. Turn off the axial fan
- 4.4.6. Wait 10-15 minutes and switch off the power to the water tank

## 5. Software operation Procedure

### 5.1. Zero Scan

- 5.1.1. Click on "ZeroScan" in the top navigation bar, then select "Auto Scan" from the drop-down menu.

### 5.2. New Method

- 5.2.1. Click on "Method" in the top navigation bar, then select "New Method" from the drop-down menu
- 5.2.2. Enter the method name and click "OK"
- 5.2.3. Set up "PARAM" as shown below, then click "OK"
- 5.2.4. Click the + icon to add an element
- 5.2.5. Enter the element name in the search box on the right side of the pop-up window, double-click the element to add its spectral line, and click 'Cancel' to close the window
- 5.2.6. Select a spectral line and set 'Repeat Standard'. Assuming there are three standard solutions, modify STDLOW to STD1 (the lowest standard) and STDHIGH to STD3 (the highest standard).
- 5.2.7. Select a spectral line and set 'Std'. Assuming there are three standard solutions, arrange the standard names in ascending order as STD1, STD2, and STD3, and input the theoretical concentration corresponding to each standard.
- 5.2.8. Select a spectral line to set 'Ms.Unit'. Click the drop-down box next to 'Ms.Unit' on the right, and choose the unit of measurement (if the required unit is not available in the drop-down box, you can manually input it).
- 5.2.9. Click 'OK' to complete the setup.

### 5.3. Select Method

- 5.3.1. Click on "**Method**" in the top navigation bar, then select "**Select Method**" from the drop-down menu.
- 5.3.2. Double-click on the method name to select.

### 5.4. Auto Search

- 5.4.1. Click on the "**Auto Search**" button at the top, or click on "**Analysis**" in the top navigation bar, then select "**Auto Search**" from the drop-down list.
- 5.4.2. Click "OK" to start.

### 5.5. Auto Attenuate

- 5.5.1. Click on the "**Auto Attenuate**" button at the top, or click on "**Analysis**" in the top navigation bar, then select "**Auto Attenuate**" from the drop-down list.

5.5.2. Click "OK" to start.

## 5.6. Standard

5.6.1. Click on the **"Standard"** button at the top, or click on **"Analysis"** in the top navigation bar, then select **"Calibration"** from the drop-down list.

5.6.2. Ensure that the standard to be tested matches the actual standard, then click the **"Measurement"** button (repeat this operation until tests for STD1, STD2, STD3... STDn are completed).

5.6.3. After completing the standard test, click the **"Curve Smoothing"** button to verify whether the working curve meets the requirements (correlation coefficient  $\geq 0.9995$ ).

## 5.7. Sample

5.7.1. Click on the **"Sample"** button at the top, or click on **"Analysis"** in the top navigation bar, then select **"Analysis Sample"** from the drop-down list.

5.7.2. Enter the sample number and click the **"Measure"** button (repeat this operation until all samples have been tested).

## 5.8. View Test Results

5.8.1. Click on **"Preparation"** in the top navigation bar, then select **"Result"** from the drop-down list

## 6. Common Fault Resolution

Process	Phenomenon	Inspection	Cause	Solution
Startup Stage	Water tank makes a beeping sound, water tank water - lack indicator light is on	Check if the water tank is short of water	Water tank is short of water	Add pure water to the water tank
	Open the right - side power supply of the instrument, the instrument does not start	Instrument main power not turned on	Voltage protection	Turn on the main power
		Gas cylinder not turned on	Pressure protection	Turn on the gas cylinder
		Gas pressure not in 0.25 - 0.4 Mpa	Pressure protection	Adjust the gas cylinder flow to 0.25 - 0.4 Mpa
		Water tank not turned on	Water pressure protection	Turn on the water tank
		Use a multimeter to check whether the water pressure switch is disconnected when the water tank is open and connected when the water tank is closed. Check if it is connected when the water tank is open and disconnected when the water tank is closed	Abnormal water pressure	Adjust the water pressure switch

Ignition Stage	no Ignition	Check the purity of argon gas	Impure argon gas	Replace with new argon gas (≥99.99)
		Ignition wire hook not hung well	Ignition wire not placed well	Stick and fasten the ignition wire hook tightly to the torch wall
		Carrier gas opened during ignition	Carrier gas will blow out the flame during ignition	Close the carrier gas before ignition
		After the high voltage is turned on, the voltage is too low	Preheating time is too short	Turn on the right - side power supply of the instrument, wait for 3 minutes of preheating, and then perform the ignition operation
		Ignition - time, there is a mesh - like upward current at the torch tube, but no flame forms	Insufficient power	After a mesh - like current appears at the torch tube, increase the voltage until the flame forms
			Mismatch	1. Open the upper right panel at the back of the instrument, find the white vacuum capacitor and make a mark
				2. Rotate the mark 1/4 turn to the left or right, and try to ignite
				3. Observe the strength of the current at the torch tube during ignition. If it becomes stronger, continue to rotate in the same direction until the flame forms; if the current becomes weaker during rotation, adjust in the opposite direction until the flame forms
				4. If the current is abnormal, adjust until the flame forms
	Ignition - after, flame goes out	Capillary tube not inserted into the solution	Air is sucked in	Ensure that the capillary tube is inserted into the solution in time
		Torch tube installed too high or too low	Improper installation position of the torch tube	Re - adjust the torch tube. Ensure that the ignition coil is between the middle layer and the outer wall of the torch tube
	Ignition - time, there is arcing between the ignition coil and the torch tube	Small water droplets on the outer wall of the ignition coil	Water droplets conduct electricity and cause arcing	Wipe off the water droplets on the coil
		The black insulating part of the ignition coil is damaged	Insulation sleeve is damaged	Replace the ignition coil
				Clean or replace the torch tube
		Torch tube touches the ignition coil	Improper matching between the torch tube and the ignition coil	Adjust the torch tube or the ignition coil to ensure that after installation, the torch tube is in the ignition coil
	Unstable Flame	Voltage and current values are too low	Insufficient power	Increase the instrument power
		Torch tube/chamber is very dirty	Solution is adsorbed on the torch tube/chamber, causing sample contamination	Soak the torch tube/chamber in a 1:1 nitric acid or hydrochloric acid solution for 24 hours, then rinse with pure water, heat to boiling for about 20 minutes, then clean with pure water and dry
				Replace the torch tube/chamber

Detection Stage		Water accumulates at the bottom of the torch tube	The interval between ignition times is too long	Dry the torch tube
		Waste liquid bottle has air bubbles	Torch tube center channel is blocked	Soak the torch tube in a 1:1 nitric acid or hydrochloric acid solution for 24 hours, then rinse with pure water, heat to boiling for about 20 minutes, then clean with pure water and dry
			Torch tube center channel is burned out	Replace the torch tube
	No such communication port or the communication port has been opened by itself	Check if the serial port is connected	Serial port line is loose	Pull out the serial port line on the computer host and reinsert it
				Pull out the serial port line on the right side of the ICP main unit and reinsert it
	Automatic alignment fails to find the peak	Capillary tube not inserted into the high - standard	Sampling error	Insert the capillary tube into the high - standard and re - align
		Capillary tube is not smooth	Nebulizer is blocked	Remove the nebulizer, open the gas path, and use fingers to blow back and forth. If the blockage still exists after multiple blows, soak the nebulizer in a 1:1 nitric acid or hydrochloric acid solution for 24 hours, then rinse with pure water and dry
		Waste liquid bottle has air bubbles	Torch tube center channel is blocked	Soak the torch tube in a 1:1 nitric acid or hydrochloric acid solution for 24 hours, then rinse with pure water, heat to boiling for about 20 minutes, then clean with pure water and dry
			Torch tube center channel is burned out	Replace the torch tube
		Room temperature changes greatly	Large temperature difference	1. Select the detection method to be used 2. Click "Modify Method" - "Analysis Line" - all lines "Restore Actual Length" 3. Re - align
				Re - configure the standard solution
	Curve fitting <0.9995	Check if the standard solution is contaminated	Standard solution is contaminated	The pure water for cleaning the capillary tube needs to be replaced daily. Note: Do not directly insert the capillary tube from the high - standard to the low - standard Note: Do not use standard solutions that have been placed for too long
		Software element concentration input error	Actual solution concentration does not match the software setting	Adjust the standard solution concentration or the concentration content in the software



		Capillary tube inlet is not smooth	Nebulizer is blocked	Remove the nebulizer, turn on the carrier gas, and block the nebulizer nozzle with your hand to blow the impurities back into the solution. If it is still blocked after back - blowing, place the nebulizer in a 1:1 nitric acid or hydrochloric acid solution and soak for 24 hours
		Small window runs peak	Instrument temperature changes too much	Control the instrument temperature to be stable at 25°C. Turn on the instrument auxiliary heating switch if necessary
		Transmission system lacks oil or is contaminated	Transmission system needs maintenance	1. Power off the instrument 2. While rotating the black knob, use a brush dipped in absolute ethanol to clean the oil stains between the copper gears twice 3. Evenly drop the watch oil between the copper gears
	Sample Result Deviation	Is the sample contaminated?	Sample is contaminated	Re - sample
		Is the sample dissolved?	Sample has precipitation	Re - sample to ensure complete dissolution
		Instrument temperature	Instrument temperature changes too much	Control the instrument temperature to be stable at 25°C. Turn on the instrument auxiliary heating switch if necessary
		Capillary tube inlet is not smooth	Nebulizer is blocked	Remove the nebulizer, turn on the carrier gas, and block the nebulizer nozzle with your hand to blow the impurities back into the solution. If it is still blocked after back - blowing, place the nebulizer in a 1:1 nitric acid or hydrochloric acid solution and soak for 24 hours
		Small window runs peak	Instrument temperature changes too much	Control the instrument temperature to be stable at 25°C.
		Transmission system lacks oil or is contaminated	Transmission system needs maintenance	1. Power off the instrument 2. While rotating the black knob, use a brush dipped in absolute ethanol to clean the oil stains between the copper gears twice 3. Evenly drop the watch oil between the copper gears
		Sample content has a large deviation from the standard curve	Sample does not match the standard	Adjust the sample concentration or re - configure the standard solution that matches the sample